Mapping Microplastics

Home About the Project How it Works Tips and FAQ Data

A citizen science initiative to map the presence of microplastic pollution in our local ecosystems.



- Baseline assessment across a wide geographic area
- Repeat testing collects information on seasonality/patterns
- Provides evidence to justify more precise/sophisticated testing

- Novel hands-on STEM activity
- Accessible to a wide audience
- Facilitates place-based learning about global plastics crisis
- Creates a pathway for local nature experiences

Evidence for using Nile Red staining to identify MP's in environmental samples:





Maes, T., Jessop, R., Wellner, N., Haupt, K., & Mayes, A. G. (2017). A rapid-screening approach to detect and quantify microplastics based on fluorescent tagging with Nile Red. Scientific Reports, 7(1), 1-10. Stanton, T., Johnson, M., Nathanail, P., Gomes, R. L., Needham, T., & Burson, A. (2019). Exploring the efficacy of nile red in microplastic quantification: A costaining approach. Environmenta I Science & Technology Letters, 6(10), 606-611. Hengstmann, E., & Fischer, E. K. (2019). Nile red staining in microplastic analysis proposal for a reliable and fast identification approach for large microplastics. Environme ntal monitoring and assessment, 191(10), 612.

Tamminga, M., Hengstmann, E., & Fischer, E. K. (2017). Nile Red Staining as a Subsidiary Method for Microplastic Quantification: A Comparison of Three Solvents and Factors Influencing Application Reliability, SDRP JESES, 2. Erni-Cassola, G., Gibson, M. I., Thompson, R. C., & Christie-Oleza, J. A. (2017). Lost, but found with Nile red: a novel method for detecting and quantifying small microplastics (1 mm to 20 µm) in environmental samples. Environmental science & technology, 51(23), 13641-13648.

Sarah Davis


Sarah Davis

Citizen Science Procedures









COLLECT

Participants use a glass jar to scoop a 500 mL sample of surface water from any local waterway.

FILTER

Participants filter the water sample using a coffee filter to remove solid debris. They photograph the sample under blue light/orange filter to document fluorescence prior to staining.


STAIN

Participants apply 1-3 mL of Nile Red stain directly onto the debris, and allow it to permeate for 45 minutes - 1 hour.

PHOTOGRAPH


Participants photograph the sample after staining under blue light/orange filter and record fluorescing microplastics greater than 1 mm.

Sarah Davis



SHARE

Participants upload their data to the project website, where it is reviewed and added to an interactive GIS map of results.



Evidence for using Nile Red staining to identify MP's in environmental samples:






Maes, T., Jessop, R., Wellner, N., Haupt, K., & Mayes, A. G. (2017). A rapid-screening approach to detect and quantify microplastics based on fluorescent tagging with Nile Red. Scientific Reports, 7(1), 1-10. Stanton, T., Johnson, M., Nathanail, P., Gomes, R. L., Needham, T., & Burson, A. (2019). Exploring the efficacy of nile red in microplastic quantification: A costaining approach. Environmenta I Science & Technology Letters, 6(10), 606-611. Hengstmann, E., & Fischer, E. K. (2019). Nile red staining in microplastic analysis proposal for a reliable and fast identification approach for large microplastics. Environme ntal monitoring and assessment, 191(10), 612.

Tamminga, M., Hengstmann, E., & Fischer, E. K. (2017). Nile Red Staining as a Subsidiary Method for Microplastic Quantification: A Comparison of Three Solvents and Factors Influencing Application Reliability, SDRP JESES, 2. Erni-Cassola, G., Gibson, M. I., Thompson, R. C., & Christie-Oleza, J. A. (2017). Lost, but found with Nile red: a novel method for detecting and quantifying small microplastics (1 mm to 20 µm) in environmental samples. Environmental science & technology, 51(23), 13641-13648.


Sarah Davis





HAWK TANK





WWW.MAPPINGMICROPLASTICS.ORG SARAHDAVIS@URI.EDU

Acknowledgments:



Dr. Kevin Matteson (Miami University) Dr. Charlie Mazel (Nightsea) Dr. Coleen Suckling (URI) TY Dr. Andy Davies (URI)

HAWK TANK





